

Development of novel immunochromatographic tests for rapid identification of OXA-48 and KPC carbapenemases in Enterobacteriaceae

C. Borlon¹, L. Denorme¹, S. Evrard², I. Ote¹, P. Bogaerts², Y. Glupczynski² and P. Mertens¹

¹Coris Bioconcept, Science Park Crealys, Rue Jean Sonet 4a, 5032 Gembloux, BELGIUM

²National Reference Laboratory for Monitoring of Antimicrobial Resistance in Gram-Negative Bacteria, CHU Dinant Godinne, Université catholique de Louvain (UCL), Yvoir, Belgium

Introduction

OXA-48 and **KPC** are two of the major carbapenemases expressed in Enterobacteriaceae and their rapid identification presents a serious challenge for diagnostic laboratories. We developed two monoplex immunochromatographic assays: **OXA-48 K-SeT** (already commercially available) and **KPC K-SeT**, and initiated investigations for a combi assay: **OXA-48 & KPC K-SeT** detecting specifically both OXA-48 and KPC carbapenemases from bacterial colonies in 15 minutes.

KPC antibodies selection

CAPTURE ANTIBODY		DETECTION ANTIBODY (Coating)										
(Gold-Conjugate)	UF4	AA6	BA4	KA5	KB4	MF5	NA3	NE7	PF7	RF8	TA1	Signal intensity
Conjugate												
UF4	-	-	-	-	-	-	-	-	+	-	-	- Negative
AA6	-	-	-	-	-	-	-	-	+	-	-	(-) Very weak positive
BA4	-	(+)	-	-	-	-	++	(+)	++	-	-	+ Weak positive
KA5	-	-	-	(+)	-	-	+	+	+	-	-	+(-) Positive
KB4	-	-	-	-	-	-	(+)	+	(+)	-	-	++ Positive/strong positive
MF5	-	-	(+)	-	-	-	-	-	+	+	-	+++ Positive/strong positive
NA3	-	(+)	(+)	+	-	-	-	-	++	-	-	****
NE7	(+)	+	+	(+)	+	-	-	-	++	-	-	
PF7	++	(+)	(+)	+	-	-	++	+++	+	+	-	
RF8	(+)	-	-	-	-	-	-	-	+	-	-	
TA1	-	-	-	-	-	-	(+)	-	+	-	-	

=> Selection of **NE7** as detection antibody and **PF7** as capture antibody

KPC K-SeT evaluation

Blind evaluation was performed on collection and clinical bacterial isolates

		Retrospective validation		Prospective validation (ongoing)		Total
		KPC K-SeT	KPC K-SeT	KPC K-SeT	KPC K-SeT	
KPC PCR (Eazyplex)	+	12	0	21	0	33
	- (Carba non KPC)	0	20	0	68	88
	- (Non carba)	0	0	0	48	48
Total		12	20	21	116	169

100% sensitivity and 100% specificity

Using pure recombinant protein, detection level of **KPC K-SeT** was defined at: **0,625 ng/ml**.

OXA-48 K-SeT performances


Blind evaluation was performed on collection and clinical bacterial isolates

		Retrospective validation		Prospective validation (ongoing)		Total
		OXA-48 K-SeT	OXA-48 K-SeT	OXA-48 K-SeT	OXA-48 K-SeT	
OXA-48 PCR	+	28	0	92	0	120
	- (Carba non OXA-48)	0	20	0	45	65
	- (Non carba)	0	28	0	51	79
Total		28	48	92	96	264

100% sensitivity and 100% specificity

Using pure recombinant protein, detection level of **OXA-48 K-SeT** was defined at: **0,125 ng/ml**.

Combi OXA-48 & KPC K-SeT development



- The best capture/detection couple of antibodies for each target (**OXA-48** and **KPC**) was selected for further combi tests development.
- Optimization tests were performed to avoid cross reactions between carbapenemases
- Buffer optimizations were done to avoid background signals

Conclusion & Future plans

We developed a new lateral flow assay (**KPC K-SeT**) for the **rapid detection** of KPC carbapenemase, and started the set-up of a combined test for the detection of both **OXA-48** (and the variants of OXA-48 family) and KPC carbapenemases (**OXA-48 & KPC K-SeT**). The **KPC K-SeT** test was further evaluated in **prospective studies in reference centers**, confirming the **excellent sensitivity and specificity of the test**. The clinical performances of the test will be further assessed on a **larger panel** of clinical bacterial isolates when we will start clinical evaluation of our combi assay.



Acknowledgments: This work was funded by the program CWALity from the Region Wallonne, Belgium

Contacts:
Scientific: pascal.mertens@corisbio.com, celine.borlon@corisbio.com
Commercial: clientcare@corisbio.com, salah.azzi@corisbio.com

