

# HAT Sero K-Set


***In Vitro* rapid detection test of *T.b. gambiense* specific antibodies in blood or plasma**

- Easy to use
- Rapid (15 min)
- On whole blood or plasma
- No refrigerators or any other equipment required
- No power supply required



[www.corisbio.com](http://www.corisbio.com)

# HAT Sero K-SeT



Performances: field external evaluation (DR Congo) on whole blood, compared to parasitology (n=493)

		95% confidence interval
Specificity	98.6%	(97.1 to 100%)
Sensitivity	98.5%	(95.9 to 100%)
PPV	96.4%	(91.3 to 98.6%)
NPV	99.4%	(97.8 to 99.9%)
Agreement	98.7%	

For more technical information

[client.care@corisbio.com](mailto:client.care@corisbio.com)


Ordering Information

Product name	Description	Storage	Code
HAT Sero K-SeT	40 tests, buffer, heparinized capillary tubes	4°C to 30°C	K-15S2
With collection set	alcohol prep pads, blood lancets		
HAT Sero K-SeT	40 tests, buffer, heparinized capillary tubes	4°C to 30°C	K-12S2

## All diagnostics available in Coris BioConcept catalogue :

Adenovirus, Adenovirus (Respiratory), Adenovirus 40 and 41, *Clostridium difficile*, Bovine Coronavirus, *Cryptosporidium*, *Escherichia coli* O157, *Escherichia coli* K99, *Giardia*, *Helicobacter pylori*, Influenza A and B viruses, *Legionella pneumophila*, Respiratory Syncytial Virus, Rotavirus, *Trypanosoma brucei gambiense*.  
Mefloquine, Proguanil.  
*Leishmania spp.*, *Trypanosoma cruzi*, *P. aeruginosa* efflux.

Distributed by:



Science Park CREALYS  
Rue Jean Sonet 4A - 5032 Gembloux - Belgium  
Ph.: +32 (0) 81 719 910 - Fx.: +32 (0)81 719 919  
sales@corisbio.com - www.corisbio.com

MADE IN BELGIUM

Fo 7.5-35-23/02  
August 2014