

Experience the colour

CE

Pylori-Strip

Rapid diagnostic test for *in vitro* detection of *Helicobacter pylori* antigen in stool specimens


- Easy to use
- Easy to read
- Specific
- Sensitive
- Reliable
- Rapid format (10')



www.corisbio.com



Accurate, economical and fast procedure!



Performance compared to EIA (n = 93)

Specificity	97.1% (33/34)
Sensitivity	96.6% (57/59)
PPV	98.3% (57/58)
NPV	94.3% (33/35)
Accuracy	96.8% (90/93)

Bibliography

For more information:

www.corisbio.com/public/science

Or send an e-mail to: info@corisbio.com

Ordering Information

Enteric range (25 tests)

Description	Code
Rota-Strip	C-1001
Adeno-Strip	C-1002
40/41 Adeno-Strip	C-1003
Combi-Strip	C-1004
Crypto-Strip	C-1005
O157 Coli-Strip (20 tests)	C-1011
Giardia-Strip	C-1013
GastroVir-Strip	C-1016
Crypto/Giardia Duo-Strip	C-1018

Respiratory range (25 tests)

Description	Code
RSV Respi-Strip	C-1006
Adeno Respi-Strip	C-1009
Influ-A Respi-Strip	C-1010
Influ-A&B Respi-Strip	C-1012
Legionella V- Test*	C-1815

Gastric range (25 tests)

Description	Code
Pylori-Strip	C-1019

All our products are also available individually packed by 10 tests.

* Only available in 10 or 20 tests.



Science Park CREALYS
 Rue Jean Sonet 4A - 5032 Gembloux - Belgium
 Ph.: +32 (0)81 719 917 - Fx.: +32 (0)81 719 919
info@corisbio.com - www.corisbio.com
 Made in Belgium

Distributed by: