




Combi K-Set


Rapid diagnostic test for *in vitro* detection of Rotavirus & Adenovirus in stool specimens

- Safe
- Easy and ready to use
- Complete set of sampling devices available
- Accurate
- Easy to interpret



www.corisbio.com

Combi K-SeT



Performance of the kit

Compared to RT-qPCR Rotavirus detection (n=194) (95% confidence interval)

Sensitivity	91.7 %	(82.1 to 96.6 %)
Specificity	99.2 %	(99.2 to 100 %)
PPV	98.5 %	(90.9 to 99.9 %)
NPV	95.3 %	(89.6 to 98.1 %)
Agreement	96.4 %	(187/194)

Compared to qPCR Adenovirus 40/41 detection(n=190) (95% confidence interval)

	82.6 %	(60.5 to 94.3 %)
	97.0 %	(92.8 to 98.9 %)
	79.2 %	(57.3 to 92.1 %)
	97.6 %	(93.6 to 99.2 %)
	95.3 %	(181/190)


For more technical information

client.care@corisbio.com

Ordering Information

PRODUCT NAME	DESCRIPTION	STORAGE	CODE
Combi K-SeT	20 cassettes, buffer	4°C to 30°C	K-1204
Combi K-SeT with collection set	20 cassettes, 20 faecal sampling systems	4°C to 30°C	K-1504
Combi-Strip	25 tests, buffer	4°C to 30°C	C-1004

— Distributed by —



CREALYS Science Park
Rue Guillaume Fouquet 11
5032 Gembloux - BELGIUM
T +32 (0)81 719 910 - F +32 (0)81 719 919
sales@corisbio.com - www.corisbio.com

MADE IN BELGIUM

D-CUS-10-04-K/04
March 2023