

40/41 Adeno-Strip


Rapid diagnostic test for *in vitro* detection of group F Adenovirus in stool specimens

- Easy to use
- Easy to read
- Specific
- Sensitive
- Reliable
- Rapid (10')



www.corisbio.com

40/41 Adeno-Strip



Performance compared to EIA (n = 153)

		(95% confidence interval)
Sensitivity	100%	(62.9 to 100%)
Specificity	99.3%	(95.6 to 100%)
PPV	90%	(54.1 to 99.5%)
NPV	100%	(96.7 to 100%)
Accuracy	99.4%	

For more technical information

client.care@corisbio.com

Ordering Information

Product Name	Description	Storage	Code
40/41 Adeno-Strip	25 tests, buffer	4°C to 30°C	C-1003

All diagnostics available in Coris BioConcept catalogue:

Virus: 40/41 Adenovirus, Adenovirus, Influenza A&B, Respiratory Adenovirus, Respiratory Syncytial Virus, Rotavirus



Bacteria: *Clostridium difficile*, *Helicobacter pylori*, *Legionella pneumophila*, *Streptococcus* (group A)

Parasite: *Cryptosporidium*, *Giardia*, *Human African Trypanosomiasis*

Antimicrobial Resistance (RESIST): OXA-23, OXA-48, OXA-163, IMP, KPC, NDM, VIM

Drug monitoring: Mefloquine, Proguanil

Distributed by:



Science Park CREALYS
Rue Jean Sonet 4A - 5032 Gembloux - Belgium
Ph.: +32 (0) 81 719 910 - Fx.: +32 (0)81 719 919
sales@corisbio.com - www.corisbio.com

MADE IN BELGIUM

D-CUS-10-03-S/07
February 2019