




# Giardia K-Set


Rapid diagnostic test for *in vitro* detection of *Giardia* in stool specimens

- Safe
- Easy and ready to use
- Complete set of sampling devices available
- Accurate
- Easy to interpret



[www.corisbio.com](http://www.corisbio.com)

# Giardia K-SeT



Performance compared to immuno-enzymatic rapid test (n = 100)

	Microscopic analysis (n=184) (95% confidence interval)		ELISA (n=94) (95% confidence interval)	
Sensitivity	89.2 %	(73.6 to 96.5 %)	100 %	(81.5 to 100%)
Specificity	99.3 %	(95.7 to 100 %)	95.8 %	(87.5 to 98.9%)
PPV	97.1 %	(83.0 to 99.9 %)	88 %	(67.7 to 96.9%)
NPV	97.3 %	(92.9 to 99.1 %)	100 %	(93.4 to 100%)
Agreement	97.3 %	(178/184)	96.8 %	(91/94)

For more technical information

[client.care@corisbio.com](mailto:client.care@corisbio.com)

Ordering Information

Product Name	Description	Storage	Code
Giardia K-SeT	20 cassettes, buffer	4°C to 30°C	K-1213
Giardia K-SeT, with collection set	20 cassettes, 20 faecal sampling systems	4°C to 30°C	K-1513

#### All diagnostics available in Coris BioConcept catalogue:

Virus: 40/41 Adenovirus, Adenovirus, Influenza A&B, Respiratory Adenovirus, Respiratory Syncytial Virus, Rotavirus



Bacteria: *Clostridium difficile*, *Helicobacter pylori*, *Legionella pneumophila*, *Streptococcus* (group A)

Parasite: *Cryptosporidium*, *Giardia*, *Human African Trypanosomiasis*

Antimicrobial Resistance (RESIST): OXA-23, OXA-48, OXA-163, IMP, KPC, NDM, VIM

Drug monitoring: Mefloquine, Proguanil

Distributed by:



Science Park CREALYS  
Rue Jean Sonet 4A - 5032 Gembloux - Belgium  
Ph.: +32 (0) 81 719 910 - Fx.: +32 (0)81 719 919  
sales@corisbio.com - www.corisbio.com

MADE IN BELGIUM

D-CUS-10-13-K/02  
February 2020