


Crypto/Giardia Duo-Strip


Rapid diagnostic test for *in vitro* detection of *Cryptosporidium* & *Giardia lamblia* in stool specimens

- Easy to use
- Easy to read
- Specific
- Sensitive
- Reliable
- Rapid (15')
- One-side stick



www.corisbio.com

Crypto/Giardia Duo-Strip



Performance of the kit

	Compared to Microscopic analysis Cryptosporidium (n=184) (95% confidence interval)		Compared to Microscopic analysis Giardia (n=184) (95% confidence interval)	
Sensitivity	86.7%	(68.4 to 95.6%)	89.2%	(73.6 to 96.5%)
Specificity	100%	(97.0 to 100%)	99.3%	(95.7 to 100%)
PPV	100%	(84.0 to 100%)	97.1%	(83.0 to 99.9%)
NPV	97.5%	(93.2 to 99.2%)	97.3%	(92.9 to 99.1%)
Agreement	97.8%	(180/184)	97.3%	(179/184)

For more technical information

client.care@corisbio.com


Ordering Information

PRODUCT NAME	DESCRIPTION	STORAGE	CODE
Crypto/Giardia Duo-Strip	25 tests, buffer	4°C to 30°C	C-1018



CREALYS Science Park
Rue Guillaume Fouquet 11
5032 Gembloux - BELGIUM
T +32 (0)81 719 910 - F +32 (0)81 719 919
sales@corisbio.com - www.corisbio.com

— Distributed by —



MADE IN BELGIUM

D-CUS-10-18-S/06
March 2023