






Giardia-Strip


Rapid diagnostic test for *in vitro* detection of *Giardia lamblia* in stool specimens

- Easy to use
- Easy to read
- Specific
- Sensitive
- Reliable
- Rapid (15')



www.corisbio.com

Giardia-Strip



Performance of the kit

Compared to ELISA (n=94)

		(95% confidence interval)
Sensitivity	100 %	(81.5 to 100 %)
Specificity	95.8 %	(87.5 to 98.9 %)
PPV	88.0 %	(67.7 to 96.9 %)
NPV	100 %	(93.4 to 100 %)
Agreement	96.8 %	(91/94)

Compared to Microscopic analysis (n=56)


		(95% confidence interval)
Sensitivity	100 %	(71.7 to 100 %)
Specificity	95.3 %	(82.9 to 99.2 %)
PPV	86.7 %	(58.4 to 97.7 %)
NPV	100 %	(89.3 to 100 %)
Agreement	96.4 %	(54/56)

For more technical information

client.care@corisbio.com

Ordering Information

PRODUCT NAME	DESCRIPTION	STORAGE	CODE
Giardia-Strip	25 tests, buffer	4°C to 30°C	C-1013
Giardia K-SeT with collection set	20 cassettes, 20 faecal sampling systems	4°C to 30°C	K-1513



CREALYS Science Park
Rue Guillaume Fouquet 11
5032 Gembloux - BELGIUM
T +32 (0)81 719 910 - F +32 (0)81 719 919
sales@corisbio.com - www.corisbio.com

— Distributed by —



MADE IN BELGIUM

D-CUS-10-13-S/17
March 2023