


HAT Sero K-Set


***In Vitro* rapid detection test of *T.b. gambiense* specific antibodies in whole blood, serum or plasma**

- Easy to use
- Rapid (15')
- On whole blood or plasma
- No refrigerators or any other equipment required
- No power supply required



www.corisbio.com

HAT Sero K-Set



Performance of the kit

		Compared to parasitology (n=493) (95% confidence interval)
Sensitivity	98.5%	(94.2 to 99.7%)
Specificity	98.6%	(96.6 to 99.5%)
PPV	96.4%	(91.3 to 98.6%)
NPV	99.4%	(97.8 to 99.9%)
Agreement	98.6%	(486/493)


For more technical information

client.care@corisbio.com

Ordering Information

PRODUCT NAME	DESCRIPTION	STORAGE	CODE
HAT Sero K-Set	40 tests, buffer, capillary tubes	4°C to 30°C	K-12S2
HAT Sero K-Set with collection set	40 tests, buffer, capillary tubes, alcohol prep pads, blood lancets	4°C to 30°C	K-15S2

— Distributed by —



CREALYS Science Park
Rue Guillaume Fouquet 11
5032 Gembloux - BELGIUM
T +32 (0)81 719 910 - F +32 (0)81 719 919
sales@corisbio.com - www.corisbio.com

MADE IN BELGIUM

D-CUS-10-S2-K/06
June 2023