


Influ A+B K-Set


Rapid diagnostic test for *in vitro* detection of Influenza A & Influenza B Viruses in nasopharyngeal secretions

- Easy to use
- Easy to read
- Specific
- Sensitive
- Reliable
- Rapid (15')



www.corisbio.com

Influ A+B K-SeT



Performance of the kit


	Compared to Immunofluorescent method			
	Influenza to A (n=100)		Influenza to B (n=100)	
Sensitivity	100%	(65.6 to 100%)	100%	(59.8 to 100%)
Specificity	100%	(94.9 to 100%)	100%	(95.0 to 100%)
PPV	100%	(65.6 to 100%)	100%	(59.8 to 100%)
NPV	100%	(94.9 to 100%)	100%	(95.0 to 100%)
Agreement	100%	(100/100)	100%	(100/100)

For more technical information

client.care@corisbio.com

Ordering Information

PRODUCT NAME	DESCRIPTION	STORAGE	CODE
Influ A+B K-SeT	20 cassettes, buffer	4°C to 30°C	K-1212
Influ A+B K-SeT with collection set	20 cassettes, buffer, 20 swabs, 20 tube-droppers	4°C to 30°C	K-1512



CREALYS Science Park
Rue Guillaume Fouquet 11
5032 Gembloux - BELGIUM
T +32 (0)81 719 910 - F +32 (0)81 719 919
sales@corisbio.com - www.corisbio.com

— Distributed by —



D-CUS-10-12-K/06
March 2023

MADE IN BELGIUM