


Pylori-Strip


Rapid diagnostic test for *in vitro* detection of *Helicobacter pylori* antigen in stool specimens

- Easy to use
- Easy to read
- Specific
- Sensitive
- Reliable
- Rapid (15')



www.corisbio.com

Pylori-Strip



Performance of the kit

Compared to Composite reference method (n=92)

(95% confidence interval)

Sensitivity	97.8%	(87 to 99.9%)
Specificity	93.5%	(81.1 to 98.3%)
PPV	93.8%	(81.8 to 98.4%)
NPV	97.7%	(86.5 to 99.9%)
Agreement	95.7%	(88/92)


For more technical information

client.care@corisbio.com

Ordering Information

PRODUCT NAME	DESCRIPTION	STORAGE	CODE
Pylori-Strip	25 tests, buffer	4°C to 30°C	C-1019
Pylori K-SeT with collection set	20 cassettes, 20 faecal sampling systems	4°C to 30°C	K-1519

— Distributed by —



CREALYS Science Park
Rue Guillaume Fouquet 11
5032 Gembloux - BELGIUM
T +32 (0)81 719 910 - F +32 (0)81 719 919
sales@corisbio.com - www.corisbio.com

MADE IN BELGIUM

D-CUS-10-19-S/07
March 2023