





RSV K-Set


Rapid diagnostic test for *in vitro* detection of Respiratory Syncytial Virus in nasopharyngeal secretions

- Safe
- Easy and ready to use
- Complete set of sampling devices available
- Accurate
- Easy to interpret



www.corisbio.com

RSV K-SeT



Performance of the kit


	Compared to FIA (n=200)	
		(95 % confidence interval)
Sensitivity	98.3 %	(89.9 to 99.9 %)
Specificity	100 %	(96.7 to 100 %)
PPV	100 %	(92.4 to 100 %)
NPV	99.3 %	(95.5 to 99.9 %)
Agreement	99.5 %	(199/200)

For more technical information

client.care@corisbio.com

Ordering Information

PRODUCT NAME	DESCRIPTION	STORAGE	CODE
RSV K-SeT	20 cassettes, buffer	4°C to 30°C	K-1206
RSV K-SeT with collection set	20 cassettes, buffer, 20 swabs, 20 tube-droppers	4°C to 30°C	K-1506
RSV Respi-Strip	25 tests, buffer	4°C to 30°C	C-1006



CREALYS Science Park
Rue Guillaume Fouquet 11
5032 Gembloux - BELGIUM
T +32 (0)81 719 910 - F +32 (0)81 719 919
sales@corisbio.com - www.corisbio.com

— Distributed by —



D-CUS-10-06-K/05
March 2023

MADE IN BELGIUM