



Detect and identify β -lactamases in AntiMicrobial Resistance FASTER from 2 types of specimens

NEW

ReSCape



Culture reagents kit for the screening of C-PGNB carriage on **rectal swabs**

- ✓ Reduced turn-around time (TAT)
- ✓ Easy-to-use
- ✓ Better infection control management in outbreak setting
- ✓ Compatible for all the targets of our RESIST kits


NEW

RESIST-BC




Culture reagents kit for CPE detection on **positive blood cultures** in patients with bloodstream infection

- ✓ Faster medical decisions
- ✓ Easy-to-use
- ✓ Help clinicians to optimize patient management
- ✓ Compatible for all the targets of our RESIST kits




✓ Rectal swab




✓ Positive blood culture




Save precious time, compared to conventional methods !!



ReSCape



RESIST-BC



DISCOVER THE WHOLE RESIST RANGE

Product name	Targets	Description	Code
OXA-48 K-SeT	OXA-48	20 cassettes, buffer, 20 tube-droppers	K-15R1
RESIST-3 O.O.K. K-SeT	OXA-48 & OXA-163 & KPC	20 cassettes, buffer, 20 tube-droppers	K-15R4
RESIST-3 O.K.N. K-SeT	OXA-48 & KPC & NDM	20 cassettes, buffer, 20 tubes and pipettes	K-15R5
O.K.N.V.I. RESIST-5	OXA-48 & KPC & NDM & VIM & IMP	2x20 cassettes, buffer, 20 tubes and pipettes	K-15R11
RESIST ACINETO	OXA-23 & OXA-40/58 & NDM	20 cassettes, buffer, 20 tubes and pipettes	K-15R13
RESIST CTX-M	Group x (CTX-M-8, CTX-M-9, CTX-M-14,...) & Group 1 (CTX-M-1, CTX-M-3, CTX-M-14,...)	20 cassettes, buffer, 20 tubes and pipettes	K-15R14

ADDITIONAL MATERIAL TO THE RESIST RANGE

Product name	Description	Code
RESIST-BC	Reagents kit for blood culture: 1 RBCL solution, 1 MS buffer, 1 Washing buffer	S-1001
ReSCape	Culture reagents kit on rectal swab: 20 CProBE MEDIUM tubes, 2 Selective Mix vials, 2 Water vials	S-1002
RESIST penta O.K.N.V.I. Control	Pack of 5 individual vials of 1 ml (4 tests per vial): OXA-48 & KPC & NDM & VIM & IMP	P-10R11
OXA-163 Positive Control	Individual control of 1 ml (4 tests) for OXA-163	P-10R4-1
OXA-23 Positive Control	Individual control of 1 ml (4 tests) for OXA-23	P-10R7
OXA-40 Antigens	Individual antigen of 1 ml (4 tests) for OXA-40	P-10R13-1
OXA-58 Antigens	Individual antigen of 1 ml (4 tests) for OXA-58	P-10R13-2
CTX-M-14 Antigens	Individual antigen of 1 ml (4 tests) for CTX-M-14	P-10R14-1
CTX-M-15 Antigens	Individual antigen of 1 ml (4 tests) for CTX-M-15	P-10R14-2
Negative Control	Individual control of 1 ml (4 tests) for all kits	CTR-1000

NEW



CREALYS Science Park
Rue Guillaume Fouquet 11
5032 Gembloux - BELGIUM
T +32 (0)81 719 910 - F +32 (0)81 719 919
sales@corisbio.com - www.corisbio.com

— Distributed by —



MADE IN BELGIUM



D-CUS-10-RES/06
December 2024