






Rota-Strip


Rapid diagnostic test for *in vitro* detection of Rotavirus in stool specimens

- Easy to use
- Easy to read
- Specific
- Sensitive
- Reliable
- Rapid (10')



www.corisbio.com

Rota-Strip



Performance of the kit

Compared to ELISA (n=214)

(95 % confidence interval)


Sensitivity	98.1 %	(92.7 to 99.7 %)
Specificity	100 %	(95.7 to 100 %)
PPV	100 %	(95.6 to 100 %)
NPV	98.2 %	(92.9 to 99.7 %)
Accuracy	99.1 %	(212/214)

For more technical information

client.care@corisbio.com

Ordering Information

PRODUCT NAME	DESCRIPTION	STORAGE	CODE
Rota-Strip	25 tests, buffer	4°C to 30°C	C-1001



CREALYS Science Park
Rue Guillaume Fouquet 11
5032 Gembloux - BELGIUM
T +32 (0)81 719 910 - F +32 (0)81 719 919
sales@corisbio.com - www.corisbio.com

— Distributed by —



MADE IN BELGIUM

D-CUS-10-01-S/09
March 2023