

Rapid detection of carbapenemase-producing Enterobacteriaceae directly from positive blood cultures by a new immunochromatographic assay

Axel Hamprecht^{a, b, #}, Harald Seifert^{a, b} and Ahmad Saleh^{a, b}

^aInstitute for Medical Microbiology, Immunology and Hygiene, University Hospital of Cologne, Cologne, Germany; ^bDZIF (German Centre for Infection Research), partner site Cologne-Bonn; [#]presenting author



UNIKLINIK
KÖLN

Institute for Medical Microbiology,
Immunology and Hygiene
University Hospital of Cologne



Background

- Bloodstream infections caused by carbapenemase-producing Enterobacteriaceae (CPE) associated with treatment failure and increased mortality
- Detection of CPE from blood cultures (BC) by standard methods takes 20-96 hours

Objectives

- to develop and evaluate a protocol for the rapid detection of carbapenemases directly from positive BC using a new multiplex immunochromatographic test (ICT)

Methods

- blood cultures spiked with Enterobacteriaceae and incubated in a BD Bactex FX system until positive
- 170 molecularly characterized clinical isolates
 - K. pneumoniae* (N=84), *E. coli* (N=53), *E. cloacae* (N=15), *E. aerogenes* (N=6), *C. freundii* (N=5), other species (N=7)
 - 126 CPE
 - 79 OXA-48-like, 18 KPC, 29 NDM
 - 44 carbapenemase negative isolates
- blood from positive BC bottles was hemolyzed, bacteria concentrated by centrifugation and lysed
- lysate was transferred to the RESIST-3 O.K.N. ICT® (Coris BioConcept, Gembloux, Belgium), which detects OXA-48-like, KPC and NDM (Fig. 1)

Funding & acknowledgements

This study was supported by grants from the faculty of medicine, University Hospital Cologne. The ICT were supplied for evaluation free of charge by Coris BioConcept.




Figure 1. Workflow for the preparation of ICT directly from positive blood cultures; SDS: sodium dodecyl sulfate; PBS: phosphate buffered saline




Figure 2. Result of the ICT from blood cultures

Results

- Without pretreatment, results cannot be easily determined because of red background
- After pretreatment with SDS, all carbapenemases can be easily identified (Fig. 2)
- Sensitivity and specificity 100% (Table 1)
- Easy and rapid protocol, time to result 20-30 min

Contact information

Axel Hamprecht, Institute for Medical Microbiology, Immunology and Hygiene, University Hospital Cologne, Goldenfelsstraße 19-21, D-50935 Cologne/Germany; Tel: (+49) 221-478 32100 e-mail: axel.hamprecht@uk-koeln.de

Results

Carbapenemase	positive tests/ total isolates	Sensitivity	Specificity
OXA-48-like	79/79	100%	100%
OXA-48	54/54		
OXA-181	8/8		
OXA-232	7/7		
OXA-162	7/7		
others	3/3		
OXA-48-like/NDM	3/3	100%	100%
KPC	18/18	100%	100%
KPC-2	17/17		
KPC-3	1/1		
NDM	29/29	100%	100%
NDM-1	27/27		
NDM-5	1/1		
NDM-7	1/1		

Table 2. Sensitivity and specificity of the ICT from blood cultures

Conclusion

- OXA-48-like, KPC and NDM carbapenemases can be reliably detected directly from positive BC bottles with the new protocol
- more rapid than other currently available assays
- can be performed in any routine microbiology laboratory
- can help to rapidly identify patients with CPE BSI and early optimize treatment